

Table Cage Sml Med Lge

- Provides light and airy conditions which suit a wide variety of flying or crawling insects.
- Suitable for use indoors or outdoors, with excellent viewing properties.
- Easily cleanable plastic base.
- Cage frame is made of rust resistant materials and can be easily assembled and disassembled.
- Has a practical curved zip access, useful for general maintenance or photography.
- Net cover is washable and replaceable.

Contents

1 x net cover 4 x rectangle sections 4 x legs 8 x anchors 8 x 12mm screws 4 x washers

Frame Assembly


- 1. Remove the screw and anchors from one end of the legs and use them to join the shaped wire sections to form a rectangle. Firmly 'hand-tight' is good enough as over-tightening may cause damage or distortion.
- 2. Remove the plastic film from the cage base.
- 3. Remove the screw and washer from the other end of the legs, leaving the anchor.
- 4. Insert a screw through each hole in the plastic base from underneath, so the screw sits in the countersunk screw hole.
- 5. Place the washer on the screw, then screw on the anchor (with leg attached). Firmly 'hand-tight' is good enough
- 6. Repeat for all the legs.
- 7. Re-fit the completed rectangle assembly on top of the legs.

Fitting the Net Cover

- 1. Slide the cover carefully over the frame, ensuring it doesn't snag.
- 2. Ease the cover underneath the base and tension the elastic cord hem, tying with a bow. .
- Access to any side of the cage can be made by easing the cover from under the base and lifting the net locally. Alternatively the long double slider zip provides a small or large general access at varying positions.

Care of the Frame and Net

- Over time our fabrics will be affected by UV rays and sunlight. When not in use disassemble and store in a cool, dry, dark place; this will maximise the life of your product. Ensure all parts are clean and dry before packing away.
- The product may be hand washed in warm (not hot) water with a mild detergent. Ensure it is thoroughly rinsed after washing. To sterilise against disease, use sterilising tablets or liquid (for example Milton®) prepared and used as per the manufacturer's instructions. Afterwards, rinse thoroughly with clean water. Repeated sterilisation may cause some minor bleaching of black net over time. Allow to air dry. Do not tumble dry. Do not iron.
- Some chemicals including phenols, alkalis, iodine and acids are likely to damage our products. This is not an exhaustive list, so if you plan to use any chemicals with our products we recommend you check a sample with the test netting we have enclosed with your order. Chemical damage can sometimes take significant take time to develop and may not be immediately obvious to the naked eye.
- Predatory or aggressive species may be encountered which are capable of chewing, biting, tearing or otherwise damaging the textile(s) of which this product is made. Wherever possible, site to minimise or remove this risk.
- If small holes are accidentally made in the fabric, a temporary repair may be made in-situ using adhesive tape pressed onto both sides of the hole. Minor damage can be repaired later by hand darning. More significant damage will require patching.
- Whilst every care is taken to ensure the surface of our plastic-coated wire and PVC board is smooth, before use we recommend that you check the board edges, ring sections and legs. If necessary, any rough areas can easily be smoothed down with fine-grade sandpaper or an emery board.



1 x base